

# On-Line Determination of Solution Supersaturation of Fast Precipitating Systems via SR ATR FTIR Spectroscopy

K.J. Roberts (U. of Leeds) and H. Groen (Heriot-Watt U.)


Abstract No. Robe9575

Beamline(s): U2B


**Introduction:** This study aims to apply SR ATR FTIR spectroscopy to an examination of the crystallisation kinetics associated with fast precipitating reactions of organic specialty chemical products, notably pharmaceuticals.

**Methods and Materials:** The experimental system comprising an AXIOM Dipper-210 ATR immersion probe connected to the NICOLET MAGNA FTIR spectrometer at beamline U2B is linked to a 400 ml computer controlled batch reactor. This enables in-situ internal reflection measurements using SR IR radiation to monitor supersaturation during precipitation and crystallisation processes.

**Results:** The experiments carried out have been directed towards characterising a poorly understood real pharmaceutical system: an organic compound, produced by Glaxo Wellcome, in aqueous solution. Calibration measurements and bench marking tests against existing lab instrumentation were successfully carried out. As can be seen in **Figure 1**, using a global source no reliable calibration curves, relating the measured ratio of two IR peaks  $RA$  to concentration of the solute and temperature, could be established. Using SR IR light as a source produced much more quantifiable results as can be seen in **Figure 2**. The fact, that the 30°C calibration curve dropped below the 15°C one, is probably due to inadequate control in experimental procedures rather than being due to the analytical technique itself.



**Figure 1.** Calibration curves for the parameter  $RA$  for different temperatures using ATR FTIR spectroscopy with a global IR source.



**Figure 2.** Calibration curves for the parameter  $RA$  for different temperatures using ATR FTIR spectroscopy with the SR IR source.

**Conclusions:** This first data on a new 2-year programme is most encouraging, particularly the data on the improved base line linearity using SR ATR FTIR. Future work will develop this work and take forward improvement in micro-probe sensor technology.

**Acknowledgments:** Work supported by EPSRC. Pharmaceutical compound provided by Glaxo Wellcome.